

# Ballert

# ORTHOPEDIC

## ORTHOTICS AND PROSTHETICS

### BALLERT OXFORD COLLAR

#### PATIENT INFORMATION

NAME: \_\_\_\_\_

ADDRESS: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

PHONE NUMBER: \_\_\_\_\_

Male: \_\_\_\_\_ Female: \_\_\_\_\_

Pt. Height: \_\_\_\_\_ Pt. Weight: \_\_\_\_\_

Mandibular Measurement: See Fig. 1(A) page 2 \_\_\_\_\_

Circumference of Neck: See Fig. 1(B) page 2 \_\_\_\_\_

Width of Neck: See Fig. 2, page 2 \_\_\_\_\_

Does weakness result in patient's head falling forward on chest or does head fall backwards?

Forward \_\_\_\_\_ Backward \_\_\_\_\_ Both \_\_\_\_\_

---

#### COST

QTY	UNIT PRICE	TOTAL
-----	------------	-------

_____	\$125.00	_____
-------	----------	-------


SHIPPING AND HANDLING:		\$12.00
------------------------	--	---------

BALANCE DUE: \_\_\_\_\_


CREDIT CARD # \_\_\_\_\_ EXP. DATE: \_\_\_\_\_

**BALLERT ORTHOPEDIC**    **PHONE (773) 878-2445**    **FAX (773) 508-6699**





**FIG. 1:** (A) Mandibular Measurement  
*Measure from jaw to center of chin.*  
 (B) Circumference of Neck



**FIG. 2:** Width of Neck Measurement

## Purpose

The Ballert® Oxford Collar we designed is a frame type cervical orthosis. It consists of a wire frame covered with foam for comfort and a spring action mandibular piece. It has a leather adjustable strap which fits behind the neck.

The purpose of the Ballert Oxford Collar is simply to help to balance the patient's head on his shoulders and prevent the neck from flexing (chin going down on chest). It is only worn for certain tasks such as reading or watching television, or for other times when patients become fatigued.



Our goal in providing orthotic devices to our patients is to improve function without interfering with the activities of daily living.

## Ballert Locations

**2434 W Peterson Ave.  
 Chicago, IL 60659  
 (773) 878-2445**

**125E. Lake Cook Road,  
 Ste 221  
 Buffalo Grove, IL 60089  
 (847) 459-9006**

**1250 N. Mill Street, Ste. 106  
 Naperville, IL 60563  
 (630) 637-9540**

**5659 S. Cottage Grove, Unit B  
 Chicago, IL 60637  
 (773) 493-2445**

**233 E. Erie St., Ste. 200  
 Chicago, IL 60611  
 (312) 787-4400**

**1725 W. Harrison Ave.,  
 Ste. 960  
 Chicago, IL 60612  
 (312) 563-2795**

**139-141 Front Street  
 Wood Dale, IL 60191  
 (630) 694-9305**

**1200 N Westmoreland Ste 100  
 Lake Forest, IL 60045  
 (773) 878-2445**